

## AFB Clinical Specimen (Mycobacteria Smear & Culture)

<b>Test description</b>	Acid fast microscopy and culture of pulmonary and non-pulmonary clinical specimens.
<b>Test use</b>	To determine the presence or absence of mycobacteria in clinical specimens and, if isolated, to identify to the species or complex.
<b>Test department</b>	Mycobacteriology Laboratory Phone: (860) 920-6649, FAX: (860) 920-6721
<b>Methodology</b>	Smear by Auramine-rhodamine fluorescent stain. Culture using agar (Lowenstein-Jensen) and broth (BACTEC™ MGIT™) media. Identification methods: DNA probe, High Performance Liquid Chromatography (HPLC), biochemical testing and growth characteristics.
<b>Availability</b>	Daily, Monday-Friday. Negative cultures are reported after 6 weeks incubation.
<b>Specimen Requirements</b>	Clinical specimens from pulmonary and extra pulmonary sites. Sputum: minimum volume 3 mL (5-10 mL preferred) collected in a sterile 50 mL conical tube. Swab specimens are strongly discouraged.
<b>Collection Kit/Container</b>	To obtain collection kit, refer to Collection Kit Ordering Information.
<b>Collection Instructions</b>	Varies according to anatomic site. If needed, site specific collection instructions may be obtained by calling the Mycobacteriology Laboratory.
<b>Specimen Handling &amp; Transport</b>	Store at 2-8° C. Transport as soon as possible to laboratory. Transport with an ice pack coolant (preferred) or at ambient temperature. Avoid temperature extremes. <u>Exceptions:</u> Room temperature storage and transport is required for blood, bone marrow, CSF, tissue, and gastric wash/lavage samples.
<b>Unacceptable Conditions</b>	Unlabeled specimens Specimens that have leaked or containers that have broken in transit Blood collected into EDTA (purple top) or ACD (yellow top tubes) Coagulated blood
<b>Requisition Form</b>	Clinical Test Requisition (select <b>AFB Clinical Specimen</b> )
<b>Required Information</b>	Name and address of submitter (and/or Horizon profile #) Patient name or identifier, town of residence (city, state, zip), date of birth. Specimen source/type, date collected, test(s) requested Please ensure patient name on the requisition matches that on the specimen.
<b>Limitations</b>	<ul style="list-style-type: none"> <li>• A negative result does not rule out infection with mycobacteria.</li> <li>• Non-acid fast organisms present in culture may interfere with isolation and identification of mycobacteria.</li> </ul>
<b>Additional information</b>	<ul style="list-style-type: none"> <li>• A nucleic acid amplification test for the presence of <i>Mycobacterium tuberculosis</i> complex DNA is automatically done on the first patient specimen submitted for AFB smear and culture that is found to be acid fast smear positive (see MTBC NAAT).</li> <li>• Anti-tuberculosis drug susceptibility testing is performed on the initial <i>M. tuberculosis</i> complex isolate from each patient.</li> <li>• Cultures are incubated for 6 weeks before being reported as negative.</li> </ul>

Revision: 8/25/15